

Automation Framework for High Intensity Focused Ultrasound Ablation Device

Development of Hybrid Automation Framework for a US-based manufacturer

Hybrid Framework

- **Data Driven Framework:** Testing with different files containing different treatment data of patients.
- **Keyword Driven Framework:** Test data and objects being independent of the automation tool.
- **Modular Automation Framework:**
 - Test Script Modularity Framework - Enables creation of Small, Independent Scripts representing Modules and Functions of the Application under Test (AUT)
 - Test Library Architecture Framework - Enables creation of library files representing modules & functions of the Application Under Test (AUT)

